



Sample Name: CL\_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   40
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   74 (F)
Injection Date  : 11/6/2025 1:18:14 PM      Inj       :   10
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.597	BV	0.3019	4.12391	1.73122e-1	0.00322
2	6.326	VB	0.3284	9.99577	4.15905e-1	0.00781
3	8.161	BV	0.3285	1.26392e5	4841.95654	98.80586
4	8.886	VB	0.0724	357.03699	79.56704	0.27911
5	13.766	BV	0.0768	140.45184	28.84316	0.10980
6	14.172	VB	0.0987	5.67330	8.40679e-1	0.00444
7	16.061	BB	0.0803	950.56146	183.99074	0.74310
8	19.781	BB	0.0765	59.69116	12.31703	0.04666

```

Totals :                      1.27919e5  5148.10422

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   40
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    74 (F)
Injection Date  : 11/6/2025 1:18:14 PM        Inj       :    10
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

\*\*\* End of Report \*\*\*