



Sample Name: CL\_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   10
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   71 (F)
Injection Date  : 11/5/2025 7:00:21 PM      Inj       :   10
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.611	BB	0.2794	3.24185	1.50468e-1	0.00240
2	8.189	BV R	0.3373	1.33192e5	4897.91699	98.75751
3	8.917	VB E	0.2935	20.90390	8.78705e-1	0.01550
4	13.773	BV	0.0769	304.99197	62.54121	0.22614
5	14.156	VB	0.0837	51.94682	9.51163	0.03852
6	16.075	BB	0.0817	1031.40356	194.94376	0.76475
7	19.789	BB	0.0746	263.23743	56.26770	0.19518

```

Totals :                      1.34868e5  5222.21046
=====

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   10
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    71  (F)
Injection Date  : 11/5/2025 7:00:21 PM        Inj       :    10
                                           Inj Volume : 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*