



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   18
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 11:53:14 PM      Inj       :    8
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.604	BV	0.2847	3.60123	1.63703e-1	0.00279
2	6.328	VB	0.3595	3.32324	1.11573e-1	0.00257
3	8.171	BV R	0.3275	1.27308e5	4768.80176	98.55018
4	8.896	VB E	0.0800	54.31924	9.89484	0.04205
5	13.809	BV	0.0833	2.37067	4.37139e-1	0.00184
6	14.147	VB	0.0872	234.88118	43.39108	0.18182
7	16.073	BB	0.0814	1055.40051	200.49127	0.81699
8	19.794	BB	0.0745	518.99139	111.02865	0.40176

```

Totals :                      1.29181e5  5134.32001

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   18
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   72 (F)
Injection Date  : 11/5/2025 11:53:14 PM      Inj       :    8
                                           Inj Volume : 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***